

Laboratory Test Updates

(effective on October 12, 2026)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	New Tests	Additional Information	Internalized Tests	Reactivated Tests	Inactive
1	C8260702	Holotranscobalamine		X	X	X	X			X		
2	C8252306	N-Telopeptides (NTx)/Creatinine Ratio, Urine				X						
3	C8252307	N-Telopeptides (NTx), 24 h urine				X						
4	C8388501	Nickel (Ni)			X							
5	C8388503	Nickel (Ni)/Creatinine Ratio, Urine			X							
6	C8378501	Manganese (Mn)			X							
7	C8378504	Manganese (Mn)/Creatinine Ratio, Urine			X							
8	C8378500	Manganese,Urine (24 Hours)			X							
9	C8265701	Beta-N-Acetyl-Glucosaminidase (NAG)				X						

1	Holotranscobalamine Test Code : 8260702 - Laboratory Services			
	Methodology : CLIA - Performed : Tue 08:00 - Reported : Wed 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum, Frozen	Separator Tube	2mL	<0°C

2	N-Telopeptides (NTx)/Creatinine Ratio, Urine Test Code : 8252306 - Laboratory Services			
	Methodology : EIA - Performed : Thu 08:00 - Reported : 17 days			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Urine, Second Morning Void	Plastic Collection Cup	5mL	<0°C

3	N-Telopeptides (NTx), 24 h urine Test Code : 8252307 - Laboratory Services			
	Methodology : EIA - Performed : Thu 08:00 - Reported : 17 days			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	24-hour Urine	24-Hour Urine Container	5mL	<0°C

4	Nickel (Ni) Test Code : 8388501 - Laboratory Services			
	Methodology : ICP-MS - Performed : Tue 08:00 Wed 08:00 Fri 08:00 - Reported : Thu 18:00 Fri 18:00 Tue 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Plasma (Heparin)	Sodium Heparin-Royal Blue Top	3mL	2-8°C

5	Nickel (Ni)/Creatinine Ratio, Urine			
	Test Code : 8388503 - Laboratory Services			
	Methodology : ICP-MS - Performed : Tue 08:00 Wed 08:00 Fri 08:00 - Reported : Thu 18:00 Fri 18:00 Tue 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Urine, Random	Plastic Collection Cup	5mL	2-8°C

6	Manganese (Mn)			
	Test Code : 8378501 - Laboratory Services			
	Methodology : ICP-MS - Performed : Tue 08:00 Wed 08:00 Fri 08:00 - Reported : Thu 18:00 Fri 18:00 Tue 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Plasma (Heparin)	Sodium Heparin-Royal Blue Top	2mL	2-8°C

7	Manganese (Mn)/Creatinine Ratio, Urine			
	Test Code : 8378504 - Laboratory Services			
	Methodology : ICP-MS - Performed : Tue 08:00 Wed 08:00 Fri 08:00 - Reported : Thu 18:00 Fri 18:00 Tue 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Urine, Random	Plastic Collection Cup	5mL	2-8°C

8	Manganese,Urine (24 Hours)			
	Test Code : 8378500 - Laboratory Services			
	Methodology : ICP-MS - Performed : Tue 08:00 Wed 08:00 Fri 08:00 - Reported : Thu 18:00 Fri 18:00 Tue 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	24-hour Urine	24-Hour Urine Container	20mL	2-8°C

9	Beta-N-Acetyl-Glucosaminidase (NAG)			
	Test Code : 8265701 - Laboratory Services			
	Methodology : ENZYMATIC/COLORIMETRY - Performed : Tue 08:00 - Reported : 2 weeks			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Urine, Frozen	Plastic Collection Cup	5mL	<0°C