

(effective on October 05, 2026)

[illegible]

1	Factor XIII (Fibrin Stabilizing Factor)			
	Test Code : 8529001 - Laboratory Services			
	Methodology : CHROMOGENIC - Performed : Mon 08:00 Fri 08:00 - Reported : Mon 18:00 Fri 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Plasma (Citrate), Frozen	Blue top (Sodium citrate)	0,5mL	<0°C

2	Antithrombin III			
	Test Code : 8530101 - Laboratory Services			
	Methodology : CHROMOGENIC - Performed : Tue 08:00 Thu 08:00 - Reported : Thu 18:00 Mon 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Plasma (Citrate), Frozen	Blue top (Sodium citrate)	0,5mL	<0°C

3	Protein C, Activity			
	Test Code : 8530301 - Laboratory Services			
	Methodology : CHROMOGENIC - Performed : Tue 08:00 Thu 08:00 - Reported : Thu 18:00 Mon 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Plasma (Citrate), Frozen	Blue top (Sodium citrate)	0,5mL	<0°C

4	Anti-Factor Xa Activity (HMW/LMW-heparin)			
	Test Code : 8552001 - Laboratory Services			
	Methodology : CHROMOGENIC - Performed : Mon-Sat 12:00 - Reported : Mon-Sat 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Plasma (Citrate), Frozen	Blue top (Sodium citrate)	0,5mL	<0°C

5	pTau-217 (Tau Pathology Biomarker)			
	Test Code : 8684922 - Laboratory Services			
	Methodology : SIMOA (Single Molecule Array) - Performed : Mon 08:00 Thu 08:00 - Reported : 1 weeks			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Plasma (EDTA), Frozen	Lavender top (EDTA)	1mL	<0°C

6	PLAC Test (Ip-PLA2)			
	Test Code : 8029953 - Laboratory Services			
	Methodology : CLIA - Performed : Wed 08:00 - Reported : Thu 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum, Frozen	Separator Tube	0,5mL	<0°C